

Experiment 32a - Compton Scattering

```
In[1]:= << CurveFit`  
CurveFit for Mathematica v7.x thru v11.x, Version 1.96, 4/4/2018  
Caltech Sophomore Physics Labs, Pasadena, CA  
  
In[6]:= SetDirectory[NotebookDirectory[]]
```

Energy calibration of MCA

Nal detector

We construct an energy calibration function for the NaI scintillator using the full-energy peaks in the calibration data. We begin by removing noise from our calibration data.

```
In[88]:= calFiles = StringJoin["32a/", #] & /@  
          {"ba133_cal_na.mca", "co60_cal_na.mca", "cs137_cal_na.mca", "na22_cal_na.mca"};  
errFiles = StringJoin["err_", #] & /@ {"ba133_cal_na.mca",  
          "co60_cal_na.mca", "cs137_cal_na.mca", "na22_cal_na.mca"};  
errFiles = StringJoin["32a/", #] & /@ errFiles;  
  
In[91]:= (* Add Poisson count errors to datafiles and export *)  
For[i = 1, i ≤ Length[calFiles], i++, (LoadFile[calFiles[[i]]];  
          MakePoisson[Xerrors -> True];  
          SaveFile[errFiles[[i]]])];  
  
In[92]:= (* Add poisson count errors to noise and export *)  
LoadFile["32a/noise_na.mca"];  
MakePoisson[Xerrors -> True];  
SaveFile["32a/err_noise_na.mca"];  
  
In[108]:= (* Removing noise by normalizing over runtime,  
and subtracting noise and propagating errors. *)  
runtimes = Table[  
          StringSplit[Import[calFiles[[i]]][[1]][[2]][[5]], {i, 1, Length[calFiles]}];  
runNoise = StringSplit[Import["32a/noise_na.mca"][[1]][[2]][[5]];
```

```

In[110]:= norm = For[i = 1, i ≤ Length[errFiles], i++, (file = Import[errFiles[[i]]];
    normFile = Table[{file[[j]][[1]],  $\frac{\text{file}[[j]][[2]]}{\text{ToExpression}[\text{runtimes}[[i]]]}$ ,
         $\frac{\text{file}[[j]][[3]]}{\text{ToExpression}[\text{runtimes}[[i]]]}$ , file[[j]][[4]]}, {j, 1, Length[file]};
    Export[StringJoin["32a/norm_", StringJoin[FileNameBase[calFiles[[i]]],
        ".dat"], normFile];)];

In[111]:= file = Import["32a/err_noise_na.mca"];
normFile = Table[{file[[j]][[1]],  $\frac{\text{file}[[j]][[2]]}{\text{ToExpression}[\text{runNoise}]}$ ,
     $\frac{\text{file}[[j]][[3]]}{\text{ToExpression}[\text{runNoise}]}$ , file[[j]][[4]]}, {j, 1, Length[file]};
Export[StringJoin["32a/norm_", StringJoin["noise_na", ".dat"], normFile];

In[112]:= (* Subtracting normalized noise from files *)
normList =
    StringJoin["32a/norm_", StringJoin[FileNameBase[#], ".dat"] & /@ calFiles;

In[113]:= normNoise = Import["32a/norm_noise_na.dat"];

In[114]:= For[i = 1, i ≤ Length[normList], i++, (file = Import[normList[[i]]];
    corrected = Table[{file[[j]][[1]], file[[j]][[2]] - normNoise[[j]][[2]],
        Sqrt[file[[j]][[3]]2 + normNoise[[j]][[3]]2},
        Sqrt[file[[j]][[4]]2 + normNoise[[j]][[4]]2}, {j, 1, Length[file]};
    Export[StringJoin["32a/corrected_", StringJoin[FileNameBase[calFiles[[i]]],
        ".dat"], corrected];)];

```

We now fit the full-energy peaks in the calibration data to construct our energy calibration function for the NaI detector.

Ba-133

```

In[115]:= With[{name = SystemDialogInput["FileOpen",
    {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}}],
    If[name != $Canceled,
        LoadFile[name]]
]

In[122]:= With[{x = SetXRange[LinearDataPlot[], Log -> False,
    Label -> "Set the X values for the range you wish to keep."}],
    Print[x];

    XRangeKeep[Sequence @@ x]
]

```

```

In[125]:= GaussianCFit[]
          calData = {}; (* {energy in keV, channel, sigchannel, errenergy=0} *)

In[127]:= calData = AppendTo[calData, {32, mean, sigmean, 0}];

In[134]:= Undo[]

In[137]:= With[{x = SetXRange[ LinearDataPlot[], Log -> False,
          Label -> "Set the X values for the range you wish to keep." ]},
          Print[x];

          XRangeKeep[Sequence @@ x]
]

```

```

In[138]:= GaussianCFit[]

In[139]:= calData = AppendTo[calData, {81, mean, sigmean, 0}];

In[145]:= Undo[]

In[149]:= With[{x = SetXRange[ LinearDataPlot[], Log -> False,
          Label -> "Set the X values for the range you wish to keep." ]},
          Print[x];

          XRangeKeep[Sequence @@ x]
]

```

```

In[150]:= GaussianCFit[]

In[157]:= calData = AppendTo[calData, {356, mean, sigmean, 0}];

```

Cs-137

```

In[158]:= With[ {name = SystemDialogInput["FileOpen",
          {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"} }}}},
          If[ name != $Canceled,
          LoadFile[name]]
]

In[165]:= With[{x = SetXRange[ LinearDataPlot[], Log -> False,
          Label -> "Set the X values for the range you wish to keep." ]},
          Print[x];

          XRangeKeep[Sequence @@ x]
]

In[166]:= GaussianCFit[]

```

```

In[168]:= calData = AppendTo[calData, {32, mean, sigmean, 0}];

In[172]:= Undo[]

In[173]:= With[{x = SetXRange[LinearDataPlot[], Log -> False,
    Label -> "Set the X values for the range you wish to keep." ]},
    Print[x];

    XRangeKeep[Sequence @@ x]
]

In[175]:= GaussianCFit[]

In[176]:= calData = AppendTo[calData, {661, mean, sigmean, 0}];

```

Na-22

```

In[179]:= With[ {name = SystemDialogInput["FileOpen",
    {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"} }}}},
    If[ name != $Canceled,
        LoadFile[name]]
]

In[180]:= With[{x = SetXRange[LinearDataPlot[], Log -> False,
    Label -> "Set the X values for the range you wish to keep." ]},
    Print[x];

    XRangeKeep[Sequence @@ x]
]

In[181]:= GaussianCFit[]

In[182]:= calData = AppendTo[calData, {511, mean, sigmean, 0}];

```

The Na-22 data seems very prone to noise, therefore we only use the 0.522 MeV annihilation peak for calibration (the others features are buried in noise. Although the non noise-corrected spectrum is cleaner, this is an indication that the calibration data collected for Na-22 is of poor quality so we neglect the other features in the calibration.).

Co-60

```

In[186]:= With[ {name = SystemDialogInput["FileOpen",
    {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"} }}}},
    If[ name != $Canceled,
        LoadFile[name]]
]

```

```
In[187]:= With[{x = SetXRange[LinearDataPlot[], Log -> False,
    Label -> "Set the X values for the range you wish to keep." ]},
    Print[x];

    XRangeKeep[Sequence @@ x]
]
```

```
In[188]:= GaussianCFit[]
```

```
In[189]:= calData = AppendTo[calData, {1173, mean, sigmean, 0}];
```

```
In[190]:= Undo[]
```

```
In[191]:= With[{x = SetXRange[LinearDataPlot[], Log -> False,
    Label -> "Set the X values for the range you wish to keep." ]},
    Print[x];

    XRangeKeep[Sequence @@ x]
]
```

```
In[192]:= GaussianCFit[]
```

```
In[193]:= calData = AppendTo[calData, {1332, mean, sigmean, 0}];
```

Constructing the calibration function

```
In[218]:= (* Combine 32 keV datapoints. *)
calData[[1]] = {32, Mean[{calData[[1]][[2]], calData[[4]][[2]]}],
    
$$\frac{\text{Sqrt}[\text{calData}[[1]][[3]]^2 + \text{calData}[[4]][[3]]^2]}{2}, 0\};$$


In[223]:= calData = Delete[calData, 4];

In[224]:= Export["32a/cal_data.dat", calData];
```

```

In[327]:= calTable = Join[{"Energy (keV)", "Channel Number"},
  Table[{StringReplace["%1", {"%1" → ToString[NumberForm[calData[[i]][[1]], 3]}],
    StringReplace["%1 ± %2",
      {"%1" → ToString[NumberForm[calData[[i]][[2]], 3]], "%2" →
        ToString[NumberForm[calData[[i]][[3]], 2]}]}], {i, 1, Length[calData]}}];
Grid[calTable, Alignment → Left, Spacings → {2, 1}, Frame → All,
  ItemStyle → "Text", Background → {{Gray, None}, {LightGray, None}}]

```

Out[328]=

Energy (keV)	Channel Number
32	38.4 ± 0.026
81	66.2 ± 0.047
356	250. ± 0.083
661	460. ± 0.092
511	358. ± 0.4
1173	834. ± 0.48
1332	944. ± 0.47

```

In[241]:= With[ {name = SystemDialogInput["FileOpen",
  {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}}],
  If[ name != $Canceled,
    LoadFile[name]]
]

```

```

In[243]:= SwitchXXandYY[]

```

```

In[244]:= With[{x = SetXRange[LinearDataPlot[], Log -> False,
  Label -> "Set the X values for the range you wish to keep." ]},
  Print[x];

  XRangeKeep[Sequence @@ x]
]



```

```

In[246]:= QuadraticFit[]

```

In[247]:= **LinearDifferencePlot[]**





$$y(x) = a + b x + c x^2$$

a=	b=	c=	
-33.3313	1.74233	-0.000707733	
$\sigma_a =$	$\sigma_b =$	$\sigma_c =$	$\chi^2 / (n-3) =$
0.119738	0.00289677	9.43175×10^{-6}	1666.42

In[248]:= **Undo[]**

```
In[249]:= With[{x = SetXRange[ LinearDataPlot[], Log -> False,  
    Label -> "Set the X values for the range you wish to keep." ]},  
    Print[x];  
  
    XRangeKeep[Sequence @@ x]  
  ]  
  
In[250]:= LinearFit[]
```


In[251]:= LinearDifferencePlot[]



$y(x) = a + b x$
 $a = 22.943$ $b = 1.38535$
 $\sigma_a = 0.535681$ $\sigma_b = 0.001081$ $\chi^2 / (n-2) = 129.9$

```
In[255]:= calibration[channel_] :=
  {Piecewise[{{-33.331281248295426` + 1.7423347894608627` channel[[1]] +
    -0.0007077325790811842` channel[[1]]^2, channel[[1]] < 300},
    {22.943046857142885` + 1.3853510522657388` channel[[1]], channel[[1]] ≥ 300}}],
  Piecewise[{{Sqrt[0.11973843181390834`^2 + channel[[1]]^2 0.002896772825801799`^2 +
    1.7423347894608627` channel[[2]]^2 + (-0.0007077325790811842` channel[[1]]^2)^2
    ( ( 9.431753627773405`*^-6 )^2 + ( 2 channel[[2]] )^2 ) ], channel[[1]] < 300},
    {Sqrt[0.5356805976904518`^2 + channel[[1]]^2 0.0010809983159092101`^2 +
    channel[[2]]^2 1.3853510522657388`^2], channel[[1]] ≥ 300}}]};
```

Energy as a function of scattering angle

The following work is done in the “Exp32.nb” notebook provided for analyzing scattering events.

40 degrees



```
In[*]:= With[ {name = SystemDialogInput["FileOpen",
  {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}}],
  If[ name != $Canceled,
    LoadFile[name]]
]

In[*]:= KeepXY[]

In[*]:= xy2x[]

In[*]:= GaussianCFit[]
```

In[]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 14.9556$ $c = 0.695793$



$\sigma_{y_{\max}} = 0.774352$ $\sigma_c = 0.632289$

$\mu = 132.636$ $\sigma = 30.3599$

$\sigma_{\mu} = 1.11627$ $\sigma_{\sigma} = 2.31373$ $\chi^2 / (n-4) = 1.34211$

```
In[ ]:= xData = {};  
In[ ]:= AppendTo[xData, {40, mean, sigmean}];  
In[ ]:= Undo[]  
In[ ]:= xy2y[]  
In[ ]:= GaussianCFit[]
```

```
In[ ]:= LinearDifferencePlot[FrameLabel -> {"Channel Number", "Count"}]
```



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 26.4203$ $c = 4.10615$

$\sigma_{y_{\max}} = 2.08407$ $\sigma_c = 1.73134$

$\mu = 354.386$ $\sigma = 15.4963$

$\sigma_{\mu} = 0.647671$ $\sigma_{\sigma} = 1.81144$ $\chi^2 / (n-4) = 1.10695$

```
In[*]:= yData = {};
```

```
In[*]:= AppendTo[yData, {40, mean, sigmean}];
```

60 degrees

```
In[*]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}},
  If[ name != $Canceled,
    LoadFile[name]]
]
```



```
In[*]:= KeepXY[]
```

```
In[*]:= xy2x[]
```

```
In[*]:= NRangeRemove[#] & /@ { 67, 68, 79, 90}
```

```
In[*]:= GaussianCFit[]
```

```
In[ ]:= LinearDifferencePlot[FrameLabel -> {"Channel Number", "Count"}]
```



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 9.29206$ $c = -1.14376$

$\sigma_{y_{\max}} = 3.66407$ $\sigma_c = 3.90427$

$\mu = 231.659$ $\sigma = 48.527$

$\sigma_{\mu} = 1.69465$ $\sigma_{\sigma} = 15.8216$ $\chi^2 / (n-4) = 1.05839$



```
In[8]:= AppendTo[xData, {60, mean, sigmean}];
```

```
In[9]:= Undo[]
```

```
In[10]:= xy2y[]
```

```
In[11]:= GaussianCFit[]
```


In[]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 24.2247$ $c = -2.69883$

$\sigma_{y_{\max}} = 2.41656$ $\sigma_c = 2.82003$

$\mu = 279.698$ $\sigma = 20.6341$

$\sigma_{\mu} = 0.692554$ $\sigma_{\sigma} = 2.62177$ $\chi^2 / (n-4) = 0.745444$

```
In[ ]:= AppendTo[yData, {60, mean, sigmean}];
```

80 degrees



```
In[ ]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}}],
      If[ name != $Canceled,
        LoadFile[name]]
    ]
```

```
In[ ]:= KeepXY[]
```

```
In[ ]:= xy2x[]
```

```
In[ ]:= GaussianCFit[]
```

In[]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 15.207$ $c = 1.69228$

$\sigma_{y_{\max}} = 0.761096$ $\sigma_c = 0.663175$

$\mu = 302.345$ $\sigma = 42.1452$

$\sigma_{\mu} = 1.36054$ $\sigma_{\sigma} = 3.19759$ $\chi^2 / (n-4) = 1.21439$



```
In[ ]:= AppendTo[xData, {80, mean, sigmean}];
```

```
In[ ]:= Undo[]
```

```
In[ ]:= xy2y[]
```

```
In[ ]:= GaussianCFit[]
```

In[]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$$y_{\max} = 58.7817 \quad c = 1.17445$$

$$\sigma_{y_{\max}} = 4.93753 \quad \sigma_c = 4.13823$$

$$\mu = 226.749 \quad \sigma_{\mu} = 14.4076$$

$$\sigma_{\mu} = 0.387345 \quad \sigma_{\sigma_{\mu}} = 1.56972 \quad \chi^2 / (n-4) = 1.04349$$

```
In[ ]:= AppendTo[yData, {80, mean, sigmean}];
```

100 degrees

```
In[ ]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}}],
      If[ name != $Canceled,
        LoadFile[name]]
    ]
```



```
In[ ]:= KeepXY[]
```

```
In[ ]:= xy2x[]
```

```
In[ ]:= NRangeRemove[#] & /@ {28}
```

```
In[ ]:= GaussianCFit[]
```

```
In[ ]:= LinearDifferencePlot[FrameLabel -> {"Channel Number", "Count"}]
```



$y(x) =$
 $y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$
 $y_{\max} = 3.08151$ $c = 1.60418$
 $\sigma_{y_{\max}} = 2.43344$ $\sigma_c = 0.876048$
 $\mu = 349.291$ $\sigma = 40.6876$
 $\sigma_{\mu} = 4.44023$ $\sigma_{\sigma} = 3649.19$ $\chi^2 / (n-4) = 0.988598$



```
In[ ]:= AppendTo[xData, {100, mean, sigmean}];
```

```
In[ ]:= Undo[]
```

```
In[ ]:= xy2y[]
```

```
In[ ]:= GaussianCFit[]
```


In[]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 19.7508$ $c = 1.8952$

$\sigma_{y_{\max}} = 1.28756$ $\sigma_c = 0.92032$

$\mu = 188.513$ $\sigma = 10.7196$

$\sigma_{\mu} = 0.574773$ $\sigma_{\sigma} = 1.24045$ $\chi^2 / (n-4) = 1.46893$

```
In[ ]:= AppendTo[yData, {100, mean, sigmean}];
```

120 degrees



```
In[ ]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}}],
      If[ name != $Canceled,
        LoadFile[name]]
    ]
```

```
In[ ]:= KeepXY[]
```

```
In[ ]:= xy2x[]
```

```
In[ ]:= GaussianCFit[]
```

```
In[ ]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]
```



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} =$
3.99783

$c =$

1.30416

$\sigma_{y_{\max}} =$

1.27369

$\sigma_c =$

0.584287

$\mu =$

360.741

$\sigma =$

56.8461

$\sigma_{\mu} =$

4.24652

$\sigma_{\sigma} =$

22.1567

$\chi^2 / (n-4) =$
0.880585



```
In[8]:= AppendTo[xData, {120, mean, sigmean}];
```

```
In[9]:= Undo[]
```

```
In[10]:= xy2y[]
```

```
In[11]:= GaussianCFit[]
```

In[]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 32.1203$ $c = -4.69935$

$\sigma_{y_{\max}} = 3.6568$ $\sigma_c = 4.22351$

$\mu = 159.316$ $\sigma = 17.1331$

$\sigma_{\mu} = 0.496917$ $\sigma_{\sigma} = 2.34467$ $\chi^2 / (n-4) = 0.96631$

```
In[ ]:= AppendTo[yData, {120, mean, sigmean}];
```

130 degrees

```
In[ ]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}}],
      If[ name != $Canceled,
        LoadFile[name]]
    ]
```



```
In[ ]:= KeepXY[]
```

```
In[ ]:= xy2x[]
```

```
In[ ]:= NRangeRemove[#] & /@ {13, 46}
```

```
In[ ]:= GaussianCFit[]
```

In[]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 4.26556$ $c = 0.496153$

$\sigma_{y_{\max}} = 4.28188$ $\sigma_c = 1.06997$

$\mu = 352.543$ $\sigma = 49.8594$

$\sigma_{\mu} = 3.65033$ $\sigma_{\sigma} = 40.9078$ $\chi^2 / (n-4) = 0.982867$



```
In[ ]:= AppendTo[xData, {130, mean, sigmean}];
```

```
In[ ]:= Undo[]
```

```
In[ ]:= xy2y[]
```

```
In[ ]:= GaussianCFit[]
```


In[]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]



$$y(x) = y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} =$	$c =$	
22.0973	1.79622	
$\sigma_{y_{\max}} =$	$\sigma_c =$	
1.56726	1.26034	
$\mu =$	$\sigma =$	
153.896	10.568	
$\sigma_{\mu} =$	$\sigma_{\sigma} =$	$\chi^2 / (n-4) =$
0.532607	1.29822	0.736339

```
In[ ]:= AppendTo[yData, {130, mean, sigmean}];
```

150 degrees

```
In[ ]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}}],
      If[ name != $Canceled,
        LoadFile[name]]
    ]
```



```
In[ ]:= KeepXY[]
```

```
In[ ]:= xy2x[]
```

```
In[ ]:= NRangeRemove[#] & /@ {75, 77}
```

```
In[ ]:= GaussianCFit[]
```

In[]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 14.07$ $c = -5.7903$

$\sigma_{y_{\max}} = 85.9889$ $\sigma_c = 72.9871$

$\mu = 330.726$ $\sigma = 120.379$

$\sigma_{\mu} = 3.60567$ $\sigma_{\sigma} = 274.349$ $\chi^2 / (n-4) = 1.08132$



```
In[8]:= AppendTo[xData, {150, mean, sigmean}];
```

```
In[9]:= Undo[]
```

```
In[10]:= xy2y[]
```

```
In[11]:= GaussianCFit[]
```

```
In[ ]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]
```



$y(x) =$

$$y_{\max} \exp\left(\frac{-(x-\mu)^2}{2\sigma^2}\right) + c$$

$y_{\max} = 56.1225$ $c = -3.08235$

$\sigma_{y_{\max}} = 2.57393$ $\sigma_c = 2.9097$

$\mu = 138.836$ $\sigma = 13.7666$

$\sigma_{\mu} = 0.357275$ $\sigma_{\sigma} = 0.95142$ $\chi^2 / (n-4) = 0.716333$

```
In[ ]:= AppendTo[yData, {150, mean, sigmean}];
```

Exporting Nal and Target Full-Energy Data

```
In[*]:= SetDirectory[NotebookDirectory[]];
In[*]:= Export["32a/na_full_energy.dat", yData];
In[*]:= Export["32a/target_full_energy.dat", xData];
```

We now recall the Nal and Target full-energy events into this notebook:

```
In[258]:= yData = Import["32a/na_full_energy.dat"];
In[259]:= xData = Import["32a/target_full_energy.dat"];
```

Testing the Compton Energy equation

We tabulate the scattered photon energy for the full energy scattering events detected above using the scatter plots, propagating all errors including the $\Delta\theta = 1^\circ$ and systematic uncertainties due to the calibration. For convenience, we obtain the total uncertainty by adding systematic and random errors in quadrature as defined in the calibration function above.

```
In[262]:= calEnergy =
  Table[calibration[{yData[[i]][[2]], yData[[i]][[3]]}], {i, 1, Length[yData]}];
In[266]:= scatteredEnergy = Table[{
   $\frac{\pi}{180}$  yData[[i]][[1]], calEnergy[[i]][[1]],
  calEnergy[[i]][[2]],  $\frac{\pi}{180}$  Sqrt[2]}, {i, 1, Length[yData]}] // N;
In[268]:= Export["32a/scattered_energy.dat", scatteredEnergy];
```

We tabulate the θ , E' data:

```

In[280]:= scatteredTable = Join[{"θ (rad)", "E' (keV)"}, Table[{StringReplace[
    "%1 ± %2", {"%1" → ToString[NumberForm[scatteredEnergy[[i]][[1]], 3]],
    "%2" → ToString[NumberForm[scatteredEnergy[[i]][[4]], 2]]}, StringReplace[
    "%1 ± %2", {"%1" → ToString[NumberForm[scatteredEnergy[[i]][[2]], 3]],
    "%2" → ToString[NumberForm[scatteredEnergy[[i]][[3]], 2]]}],
    {i, 1, Length[scatteredEnergy]}]];
Grid[scatteredTable, Alignment → Left, Spacings → {2, 1}, Frame → All,
    ItemStyle → "Text", Background → {{Gray, None}, {LightGray, None}}]

```

Out[281]=

θ (rad)	E' (keV)
0.698 ± 0.025	$514. \pm 1.1$
1.05 ± 0.025	$399. \pm 1.5$
1.4 ± 0.025	$325. \pm 0.98$
1.75 ± 0.025	$270. \pm 1.$
2.09 ± 0.025	$226. \pm 0.85$
2.27 ± 0.025	$218. \pm 0.88$
2.62 ± 0.025	$195. \pm 0.66$

We can now test the Compton energy relationship with scattering angle:

```

In[269]:= With[ {name = SystemDialogInput["FileOpen",
    {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}},
    If[ name != $Canceled,
        LoadFile[name]]
    ]

(* Transforming data to linear form *)

```

```


In[270]:= xnew[ x_, y_ ] := Sin[ $\frac{1}{2}x$ ]2
ynew[ x_, y_ ] :=  $\frac{1}{y}$ 

```


```
DataTransform[]
```

```
In[275]:= LinearFit[]
```

```
In[277]:= LinearDifferencePlot[FrameLabel -> {"Sin2( $\frac{1}{2}\theta$ )", " $\frac{1}{E'}$  (keV-1)"}]
```



```
Out[277]=
```



$$y(x) = a + b x$$

$$a = 0.00149672 \quad b = 0.00385771$$

$$\sigma_a = 0.0000289197 \quad \sigma_b = 0.0000458586 \quad \chi^2 / (n-2) = 1.6505$$

We obtain a good linear fit, indicating that the Compton relationship indeed holds. Furthermore,

```
In[285]:= k0 =  $\frac{1}{a} \{1, \frac{\text{sigma}}{a}\}$  (* keV *)
```

```
Out[285]= {668.126, 12.9096}
```



```
In[286]:= me =  $\frac{2}{b} \left\{1, \frac{\text{sigb}}{b}\right\} (* \text{ keV } *)$ 
```

```
Out[286]:= {518.442, 6.16299}
```

We therefore have estimates for $k_0 = 668 \pm 13$ keV and $m_e = 518.4 \pm 6.2$ keV, which are within error of the expected values for the full-energy photons for Cs-137 and the rest mass of the electron, respectively. This again indicates consistency with the Compton energy relationship with scattering angle.



Target scintillator calibration


Unfortunately, only the Cs-137 spectrum has a discernible, single Compton edge. The Na-22 and Ba-133 calibration data have noise issues making it impossible to accurately identify the Compton edge, and the Co-60 spectrum has two full-energy events, making it difficult to distinguish the superposed Compton edges.

We use the ComptonSpectra1.nb notebook from Experiment 30b and fit the Cs-137 Compton edge:

```
In[ ]:= ComptonFit[0.661, 5000, 450]
```

```
In[ ]:= LinearDifferencePlot[FrameLabel → {"Channel Number", "Count"}]
```





$y(x)$ = Compton spectrum from a photon with energy k_0 , scaled in X and Y and convolved with the scintillator resolution e_{pe} (energy/photoelectron).

$K0$ (MeV) =	$Tedge$ (MeV) =	
0.661	0.476728	
EPE (eV/photoelectron) =	$sigEPE$ =	
9706.37	360.519	
$Xscale$ (channels/MeV) =	$sigXscale$ =	
953.657	1.8013	
$Yscale$ =	$sigYscale$ =	
(counts/channel) / (probability/MeV)	0.618745	
102.133		
$Yoffset$ (counts/channel) =	$sigYoffset$ =	
2.50152	0.221614	
$Xedge$ (channel) =	$sigXedge$ =	χ^2 =
454.635	0.858731	0.970863

We obtain a reasonably good fit, and a rough estimate of the calibration given by Xscale above.

```
In[293]:= calibrationTarget[channel_] :=
  channel[[1]]
  953.6573547506756` * 1000 {1, Sqrt[ (1.8013021195385253`^2 + (channel[[2]] / 953.6573547506756`)^2 )]};

In[294]:= calEnergyTarget = Table[
  calibrationTarget[{xData[[i]][[2]], xData[[i]][[3]]}], {i, 1, Length[xData]};

In[296]:= targetEnergy = Table[{ (pi/180 xData[[i]][[1]], calEnergyTarget[[i]][[1]],
  calEnergyTarget[[i]][[2]], (pi/180 Sqrt[2])}], {i, 1, Length[xData]}] // N;

In[297]:= Export["32a/target_energy.dat", scatteredEnergy];
```

We tabulate the θ , E data:

```

In[302]:= targetTable = Join[{" $\theta$  (rad)", "E (keV)"}, Table[{StringReplace[
    "%1  $\pm$  %2", {"%1"  $\rightarrow$  ToString[NumberForm[targetEnergy[[i]][[1]], 3]],
    "%2"  $\rightarrow$  ToString[NumberForm[targetEnergy[[i]][[4]], 2]]}, StringReplace[
    "%1  $\pm$  %2", {"%1"  $\rightarrow$  ToString[NumberForm[targetEnergy[[i]][[2]], 3]],
    "%2"  $\rightarrow$  ToString[NumberForm[targetEnergy[[i]][[3]], 2]]}],
    {i, 1, Length[targetEnergy]}]];
Grid[targetTable, Alignment  $\rightarrow$  Left, Spacings  $\rightarrow$  {2, 1}, Frame  $\rightarrow$  All,
    ItemStyle  $\rightarrow$  "Text", Background  $\rightarrow$  {{Gray, None}, {LightGray, None}}]

```

Out[303]=

θ (rad)	E (keV)
0.698 ± 0.025	$139. \pm 1.2$
1.05 ± 0.025	$243. \pm 1.8$
1.4 ± 0.025	$317. \pm 1.5$
1.75 ± 0.025	$366. \pm 4.7$
2.09 ± 0.025	$378. \pm 4.5$
2.27 ± 0.025	$370. \pm 3.9$
2.62 ± 0.025	$347. \pm 3.8$

Consistency with conservation of energy

We add the scattered photon energy and the outgoing electron energy for each angle, and expect to obtain data points consistent with a constant at 661 keV if the interactions are consistent with energy conservation.

```

In[304]:= totalEnergy = Table[
    {scatteredEnergy[[i]][[1]], scatteredEnergy[[i]][[2]] + targetEnergy[[i]][[2]],
    Sqrt[scatteredEnergy[[i]][[3]]2 + targetEnergy[[i]][[3]]2],
    scatteredEnergy[[i]][[4]]}, {i, 1, Length[scatteredEnergy]}];



In[306]:= Export["32a/conservation.dat", totalEnergy];

In[307]:= With[ {name = SystemDialogInput["FileOpen",
    {DataFileName, {"data files"  $\rightarrow$  {"*.dat", "*.mca"}, "all files"  $\rightarrow$  {"*"}}}}],
    If[ name  $\neq$  $Canceled,
        LoadFile[name]]
]

In[309]:= LinearFit[]

```

```
In[317]:= LinearDifferencePlot[FrameLabel → {"Angle (rad)", "Total Energy (keV)"}]
```



```

y(x) = a + b x
a= 696.317      b= -48.6362
σa= 2.90004    σb= 1.99282    χ2 / (n-2) = 27.7605


```

This is clearly not consistent with a constant, indicating that the larger angle data points either contain significant systematic errors, or some of the larger angle data points are skewed by the presence of Compton scatters in the NaI scintillator. This would make sense, seeing as we found it much more difficult to extract the full-energy events from the scatter plots for the large-angle measurements. We remove the large-angle measurements and attempt the same fit.


```
In[318]:= With[{x = SetXRange[ LinearDataPlot[], Log -> False,  
    Label -> "Set the X values for the range you wish to keep." ]},  
    Print[x];  
  
    XRangeKeep[Sequence @@ x]  
]
```

```
In[322]:= LinearFit[]
```

```
In[323]:= LinearDifferencePlot[FrameLabel → {"Angle (rad)", "Total Energy (keV)"}]
```



```
Out[323]=
```



$$y(x) = a + b x$$

a=	b=	
662.547	-15.4419	
$\sigma_a =$	$\sigma_b =$	$\chi^2 / (n-2) =$
3.54747	3.20908	2.61749


This is still not consistent with a slope of 0 (constant) but we have a much better fit with a much smaller slope, with the intercept within error of the Cs-137 full-energy events. This indeed indicates as suspected earlier, the presence of systematic errors in the larger-angle measurements due to Compton scattering in the NaI scintillator difficult to extract from the full-energy events when converting the scatter plot to energy peak data.

Appendix - Calibration spectra and scatter plot

Typical NaI scintillator calibration spectrum - Cs-137

```
In[329]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}},
  If[ name != $Canceled,
    LoadFile[name]]
]
```


```
In[331]:= LinearDataPlot[FrameLabel -> {"Channel Number", "Count"}]
```



Typical NaI scintillator Compton scattering spectrum (80°) - Cs-137

```
In[334]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}},
  If[ name != $Canceled,
    LoadFile[name]]
]
```


```
In[335]:= LinearDataPlot[FrameLabel -> {"Channel Number", "Count"}]
```



Typical Plastic scintillator calibration spectrum - Cs-137

```
In[332]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}},
  If[ name != $Canceled,
    LoadFile[name]]
]
```

```
In[333]:= LinearDataPlot[FrameLabel -> {"Channel Number", "Count"}]
```




Typical Plastic scintillator Compton scattering spectrum (80°) - Cs-137

```

In[336]:= With[ {name = SystemDialogInput["FileOpen",
      {DataFileName, {"data files" -> {"*.dat", "*.mca"}, "all files" -> {"*"}}}],
  If[ name != $Canceled,
    LoadFile[name]]
]


In[337]:= LinearDataPlot[FrameLabel -> {"Channel Number", "Count"}]

```




Typical Compton scatter plots - Cs-137

40°



80°



150°

